

HTA® Hysteroscopic Sheath



- 7.8mm O.D.
- Insulated polycarbonate
- Accepts $\leq 3\text{mm}$ hysteroscope

HTA® Treatment



- Dilate to 8mm / 24 Fr.
- Insert Sheath under direct visualization
- Perform diagnostic hysteroscopy
- Position Sheath inside the internal cervical os
- Treatment at 90°C for 10 minutes
- One minute cool flush

*No Laminaria

Fluid Reservoir Height Above Uterus

cm	mm Hg
115	50
122	55
127	58
137	64
147	70*

*minimum typical tubal opening pressure

The HTA Procedure



HTA Procedure – Time lapse
10 minutes → 20 seconds