


Remember scope alignment !!!



Visual Guidance for Access




Direction of view

The 25° direction of view together with the 85° wide angle field of view allows full uterus visualization by merely rotating the scope.

Minimizes trauma to the cervix and reduces the risk of foosing distension.





Distention Medium

- IV Pole
- Pressure Bag
- 3 Liter Normal Saline- (Warm)
- Cystoscopy tubing

The photograph shows a metal IV pole with a pressure bag and a 3-liter normal saline bag hanging from it. The equipment is set up for use in a clinical setting.